LANGUAGES

English	2
Spanish	8
Italian	14
Control changes	20

INSTRUCTIONS FOR USE

This kit is for the collection and preservation of samples until they arrive at the Igenomix facilities. These instructions are for the embryologist. The kit will transport embryo samples. The materials included in this kit can only be used by medical professionals, always in medical settings. This is not a self-assessment kit. This kit is not automated.

The kit will be used for the following diagnostic tests: **PGT-SR, Smart PGT-A, Smart PGT-A/SR Plus and PGT-M.**

Basic information

- **PGT-SR (Preimplantation Genetic Testing for Structural Rearrangements)**

A genetic analysis appointment with Igenomix (labpgsspain@igenomix.com) is required first to assess the suitability and method of a PGT-SR. The possible outcomes are:

- Instant PGT-SR, with no need for prior study.
- The need for a prior study to assess the suitability of the test (pre-PGT-SR). In this case, we will need a blood sample from the carrier in an EDTA tube. We will perform NGS on this sample and see if it is possible to detect potential imbalances in the embryos. (EDTA tube not included in this kit).

- **Smart PGT-A (Preimplantation Genetic Testing for Aneuploidies)**

- **Smart PGT-A/SR Plus (Preimplantation Genetic Testing for Aneuploidies, Structural Chromosomal Rearrangements, Ploidies and embryo quality control)**

The Smart PGT-A/SR Plus is an optional analysis that adds the evaluation of ploidy and embryo quality control. The triploid and haploid differences are the most common ploidy abnormalities and are incompatible with normal growth and development. Embryo quality control includes the analysis of the relatedness between samples in a cohort and the presence of external DNA contamination.

- **PGT-M (Preimplantation Genetic Testing for Monogenic Disorders)**

A genetic analysis appointment with Igenomix (pgdmonogenicspain@igenomix.com) is required to assess the suitability and method of a PGT-M

- The need for a prior study (pre-PGT-M). The results of the genetic analysis will indicate the blood/saliva/buccal swab samples required for performing the test. Igenomix will provide a kit for the purpose of obtaining these samples.

Call +34 963 905 310 or contact us by email.

Smart PGT-A, PGT-SR and Smart PGT-A/SR Plus: supportspain@igenomix.com, labpgsspain@igenomix.com or embryology@igenomix.com

PGT-M: supportspain@igenomix.com, pgdmonogenicspain@igenomix.com or embryology@igenomix.com

- Request more information about the test
- Resolve any queries about the report
- Acquire the test
- Request printed versions or information in other languages at no extra cost.

- Before beginning a PGT cycle, it is important to have the correct documents prepared. This depends on the type of service requested:

- No biopsy or tubing services required:
 - For fresh PGT transfers: Complete the document (request form) before the embryo biopsy and send it by email to **labsamples@igenomix.com, labpgsspain@igenomix.com and embryology@igenomix.com.**
 - For deferred PGT transfers: you are not required to send the request form by email prior to the biopsy.

- Biopsy service required: send the request form with as much information as possible when you are about to start the patient COS (Controlled Ovarian Stimulation).

If you do not have the request form, please send the following information as soon as possible to **embryology@igenomix.com** (at least 1 week before the biopsy):

- Unique patient identifier/Clinical record number + patient initials.
- Start date of ovarian stimulation.
- Possible ovarian puncture date.

Include all the documents (request form + biopsy form) inside the box provided by Igenomix.

INSTRUCTIONS FOR USE

Step 2: Receiving the KIT

Check the contents of the PGT kit (it is VERY IMPORTANT to handle the contents with gloves to prevent contamination) and ensure it contains:

- 1 small box containing all the kit items:
 - 3 tubes with washing and tubing solution (washing/loading).
 - 20 PCR tubes.

The kit can arrive at room temperature. Check the expiry date on the tubes and then store the box in a fridge at 4 °C until you are going to use it.

Note: If the expiry date has elapsed or the package is damaged, dispose of it in accordance with current local regulations and contact Igenomix to get a new kit.

Step 3: Sample collection

Warnings/Precautions:

The following may lead to our laboratory rejecting the sample:

- Documents (request and informed consent form) are not properly completed.
- The mandatory fields marked with an asterisk (*) have not been filled in.
- The signature of the patient and the requesting doctor is missing on the request and informed consent form.
- The karyotype is not included on the request form in the case of PGT-SR requests.
- Samples are incorrectly labelled.
- Use only the washing and tubing solution provided in the Igenomix kit. Using other solutions may affect the accuracy of the test and impact the results.

The thermally insulated box which comes with this kit should only be used for transporting samples.

Please include the identification details of the person collecting the sample.

PGT-SR: Only admits biopsy of the blastocyst on day 5/6/7 of development and deferred transfer: biopsy of 4-8 trophectoderm cells required. The biopsy can only be made on day 3* and a fresh transfer in the case of PGT-SR for Robertsonian translocations.

Smart PGT-A: Admits both single-cell biopsy on day 3* of development and trophectoderm biopsy.

*Check suitability with your local Igenomix laboratory.

Smart PGT-A/SR Plus: Only admits biopsy of the blastocyst on day 5/6/7 of development and deferred transfer.

PGT-M: Admits biopsy of the blastocyst on day 5/6/7 of development: biopsy of 4-8 trophectoderm cells. Ask about suitability of biopsy on day 3.

See tubing procedure in Annex 1.

Note: Igenomix recommends contacting authorized waste disposal companies for the correct disposal of medical products such as leftover solutions (washing/tubing) and empty containers (washing/tubing solution containers). Please ensure to comply with current local regulations in your place of work.

Step 4: Transporting the sample

Ensure you have all the items Igenomix provides for transport:

- 1 thermally insulated box
- 1 padded envelope
- 2-3 red ice packs (store at -20 °C, at least 24 hours before sending the sample)
- 1 cool rack inside the envelope (store at -20 °C). Do not remove the rack from the envelope until ready to use.

INSTRUCTIONS FOR USE

To prevent the sample stability from being compromised, which could lead to a test failure, we recommend following these instructions:

1. Identify the samples as explained in the Tubing Procedure (Annex 1).
2. Insert the cool rack with all the samples inside the freezer (-20 °C) for at least one hour before shipping.
3. If the samples are to be shipped on a different day to the biopsy sampling, we recommend keeping the rack and samples at -20 °C until shipping.
4. Insert the previously frozen rack with the samples inside the white padded envelope provided by Igenomix. Place the frozen rack between the ice packs inside the thermally insulated box, in the same way it was sent by Igenomix. To keep the rack frozen during shipping, it is very important to use the ice packs provided by Igenomix and not any others.
5. Fold and place the signed and completed informed consent and request forms inside the box for shipping.

Step 5: Sample pick-up

For picking up the sample in Spain, the EU and other countries:

 **+34 965 023 460**  **pickup@igenomix.com**

To prevent the sample stability from being compromised, which could lead to a test failure, we recommend following these instructions:

Pick-ups in Spain:

- Monday to Friday
- A 3-hour time margin is required to arrange for sample pick-ups.

Pick-up in the EU and other countries:

- Check availability with our customer service department.
- A 24-hour time margin is required to arrange for sample pick-ups.

For picking up the sample in Italy:

 **+39 0424472449**  **supportitalia@igenomix.com**

Collection in Italy:

- Monday to Friday.
- A 24-hour time margin is required to arrange for sample pick-ups.

Holidays: The aforementioned days may differ in the event of holidays. In such cases, our customer service department will notify you in advance of the available pick-up days.

We recommend sending samples as quickly as possible in order to ensure the quality of the sample.

Step 6: Test results

Igenomix will send the test results to the email address or addresses the clinic has written on the request form.

• PGT-A and PGT-SR samples with deferred transfer:

Test results will be available 7 business days after Igenomix receive the sample.

• PGT-A and PGT-SR samples with fresh transfer (only Robertsonian translocations):

Test results will be available the morning after the day Igenomix receives the samples.

• PGT-M samples:

Test results will be available 10 business days after Igenomix receive the sample.

• Smart PGT-A Plus and PGT-SR Plus:

Test results will be available in a time depending on the type of study requested.

INSTRUCTIONS FOR USE

Additional information

In the event of serious incidents:

Any serious incident that may occur in relation to the materials included in the kit supplied must be reported to **vigilance@igenomix.com** and to the competent health authority in the Member State where the user and/or patient is located.

Symbols key:

Manufacturer

Date of expiry

European Conformity

Medical device for *in vitro* testing

Keep away from sunlight

Temperature limit

Do not use if the package is damaged

Contains biological material of human origin

Do not reuse

Keep dry

Check the instructions for use.

Annex 1: TUBING PROCEDURE

For PGT (blastomere or trophectoderm)

Equipment, instruments and reagents required (for PGT-A, PGT-SR and PGT-M):

- Laminar flow cabin (preferably with UV light).
- Magnifying glass for cell manipulation.
- Automatic pipette: 1µL-10µL. We recommend using one exclusively for this PGT procedure to prevent contamination. Otherwise, please sterilize with decontamination solution.
- Filtered sterile tips, where possible DNase/RNase-free.
- Micromanipulation pipette, such as Stripper or Flexipet to perform manual tubing with volume measurement.
- Polycarbonate capillary pipette tubes of diameter between 100-150µm [such as: VITROMED Denuding Pipettes - V-DEN or RI (Origio) - EZ-Tip®].
- 60mm-Petri dishes or other dishes suitable for ICSI
- 0.2ml PCR tubes – provided by IGENOMIX.
- Cool rack – provided by IGENOMIX.
- COLD mineral oil.
- DNA AWAY™ Surface decontamination solution (or alternative such as: Oosafe, Fermacidal, Soap 7X, etc.).
- Washing/tubing solution (Washing/Loading solution) – provided by IGENOMIX.

Method:

To avoid contamination:

- We recommend the tubing procedure be performed in a laboratory separate from, but close to, the IVF laboratory, which has its own ventilation system. In this case, all cleaning procedures must be carried out with 70% ethanol solution as decontamination solution.
- If it were not the case, non-embryo toxic decontamination solutions should be used for cleaning the material.
- The embryologist must wear personal protective equipment (mask, gloves, hair/beard cover). As little skin as possible must be exposed in the hood to prevent DNA contamination.
- Disinfected or sterile gloves must always be used when handling materials and samples.
- We recommend using equipment and materials used exclusively for PGT.
- Disposable equipment must always be freshly opened or specially reserved for PGT (always handled with gloves).
- We recommend handling the samples in clean rooms, where possible separate from other activities.
- Avoid processing DNA in this room, as well as in the cooler/freezer where the buffers are stored (washing solutions).
- Clean all the material well with decontamination solution (always using gloves) taking extra care to prevent DNA contamination.

Preparation:


Before beginning the process, we recommend you have all the material prepared so tubing can be performed immediately after performing the embryo biopsy.

- All non-disposable material used in the hood tubing area must be cleaned (always using gloves) with the chosen decontamination solution – preferably DNA AWAY™, or Oosafe/Fermacidal – as well as the hood surface and magnifying glass.
- All disposable equipment must be freshly opened or for exclusive PGT use so long as it has always been handled with gloves.
- If UV light is available, all material (NOT THE EMBRYOS) must remain exposed to the light for 5 minutes before the procedure.
- Prepare the PCR tubes depending on the number of embryo biopsies, filling them with 2.5µL of the Washing/Loading buffer provided by Igenomix (we recommend you prepare additional tubes).
- Shake the PCR tubes with the Washing/Loading buffer, checking there are no bubbles in the solution. Ensure the buffer volume remains at the bottom of the tube. If no centrifugation methods are available, vigorously shake the tubes manually.
- If you are not able to perform the tubing immediately, place the PCR tubes in a cool rack and keep them in the freezer (-20°C) until the tubing process begins. If tubing is performed immediately, place the tubes in the hood at room temperature.
- DO NOT label the tubes until the biopsy is finished, but do so before starting with the tubing process.

Annex 1: TUBING PROCEDURE

For PGT (blastomere or trophectoderm)

- Washing plate: prepare a 60mm Petri dish (or an ICSI one) according to the diagram below. You will need a row of 4 drops for each biopsy sample. Each drop must contain 10µL of the washing buffer (Washing/Loading) provided by Igenomix, which is then covered with cold mineral oil. There can be a maximum of 6 rows per dish (corresponding to 6 samples).

WASHING PLATE

Sample collection:

1. Perform the procedure at room temperature, the hood surface must be cold.
2. After the biopsy and immediately before tubing, the 0.2ml PCR tubes must be labelled using a marker with the patient initials and embryo number.
3. If the 0.2ml PCR tubes are already prepared and stored in the freezer, take them out and put them in the hood so they thaw just before tubing. Ensure that the loading buffer is in liquid form and there are no bubbles.
4. Apart from the washing plate, prepare another plate at the time of tubing with a large drop (~100µL) of washing media but do not cover with oil. This drop will be used to clean the capillary pipette tube each time you place a sample. Therefore, you should prepare 1 for each 1-2 embryo biopsies.
5. Use a new capillary pipette tube for each individual cell (day 3 biopsy) or per trophectoderm biopsy.
6. Each blastomere or trophectoderm sample should be rinsed in the first 3 drops of each row of the washing plate. Load the capillary pipette tube with clean washing buffer from the large drop each time you aspirate the sample. Place the sample in the 10µL washing buffer drop and discard the remaining solution in the capillary pipette before loading it once again with clean solution.
7. When you get to the third drop, repeat the clean buffer loading process in the capillary pipette, aspirate the sample and insert it into the PCR tube, loading it under the magnifying glass in the 2.5µL deposited at the bottom of the tube. The embryologist should see how the biopsy enters the tube under the magnifying glass. This is a key stage in the process and the biopsy must be released inside the tube with as little washing buffer as possible (maximum 1.5µL).
8. The fourth drop of each row is reserved for checking the capillary pipette in the event that the embryologist does not see the sample being released into the PCR tube. If the biopsy appears in the fourth drop, discard the PCR tube where you were supposed to have inserted the biopsy and load the sample in a newly-labelled PCR tube.
9. As you collect the samples, place them in the cool rack and complete the biopsy form matching the identification details to the samples in the PCR tubes. This form must be sent with the samples to Igenomix.
10. For PGT-M cases, we request a blank per embryo for day-3 biopsies. For trophectoderm samples, it is enough to perform one blank per case and day.
11. For PGT-A, we recommend sending one blank per case.
12. Creating blanks: after releasing each sample into the PCR tube and checking it does not remain in the capillary pipette, aspirate a small amount of the washing media from the last drop where the biopsy was released and insert it into a new tube.
13. Close the tubes well and place them in the rack.
14. Insert the cool rack with all the samples inside the freezer (-20°C) for at least one hour before shipping. If the samples are to be shipped on a different day to the biopsy sampling, we recommend keeping the rack containing the samples at -20°C until shipping.
15. To prepare for shipping, place the rack in the padded envelope and place it in the thermally-insulated box between the 2 ice packs, just as it came. Do not forget to include the PGT form in the box.

Note: Igenomix recommends contacting authorized waste disposal companies for the correct disposal of medical products such as leftover solutions (washing/loading) and empty containers (washing/loading solution containers). Please ensure to comply with current local regulations in your place of work.

If you have any questions, please contact embryology@igenomix.com.

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